


LVL'D UP

The next Generation



Redesigned 2D tubes

- 1 2D bottom code / Protective rim against mechanical influences (scratches)
- 2 1D side code / HR side code
- 3 Open space for customized coding
- 4 SAFE® low leachables medical polypropylene



New cap design

Excellent sealing
Easy manual handling
Automation-friendly cap

Your free sample pack

Email us at info@autogen.com




Experience our new LVL Racks



Easy Lock

One simple move to open and close



Anti-freeze gradient

Prevents ice formation



Tube Centering

Automation-friendly beveled domes



Curved well rim

Safety guides the tube and reduces collisions